


AAC Kit Inventory of Equipment

West Dunbartonshire



- **SB 10 Mounted and stand 000065 mains adaptor and lead, manual and pocket instructions**
- **A tabletop stand is also available from the kit for this device**
- **More information from www.sensorysoftware.com**
- **Grid 2 is the communication software package on this device**

NOVA Chat 10



NOVA chat is a communication system that offers Chat software on a 10" android platform.

Size of grid varies from a very simple 4 grid to more complex communication grids. A combination of words and symbols are used on the larger grids.

Clients should generally have a good knowledge of symbols and/or written words, be able to categorise and navigate around a low tech system before a high tech system like this would be considered.

A quick guide for using the NOVA Chat is available from

<http://cache.salttillo.com/downloads/nova-quick-reference-guide-final.pdf>

I Pad




- I Pads with apps Proloquo2go, P2GText, Pictello, Touch Chat, Go Talk Now Plus, Verbally, Speak it and Predictable
- Range of apps using symbols and/or written word.
- Client will need to be able to “scroll” “swipe” “pinch” for some of the apps
- Support for use of the apps is available through online guides

Lightwriter SL40

www.toby-churchill.com



- 000091 A, main component, B carrying case, C,D keyguard, F,E single switches, 000091H special screwdriver and mains adaptor
- User guide, SMS and Sim card



- www.attainmentcompany.com
- Serial number 0000100
- Recordable messages
- Device can take 5 levels of messages
- Top 3 messages are the same on each level and should be things that are needed the most e.g. toilet, help, hello
- Grids for the GoTalk can be made on Boardmaker
- This is recorded speech so all the messages to be accessed will need to be recorded prior to giving it to the client.

Go Talk 20+

- www.AttainmentCompany.com
- 0000108
- As for Go Talk 9 but with more messages



Go Talk 32

www.AttainmentCompany.com

000072 Includes adaptor, strap and 2legs to hold up at an angle



Step Pad



This organises recorded messages for auditory prompts or expressive communication. It helps people who have difficulty completing multistep activities/tasks without support.

More information and a user guide can be found on
www.AttainmentCompany.com

Big Step by Step



- RTS 000126
- www.ablenetinc.com
- 4 different coloured covers and lead
- Recordable device
- Can take a series of messages
- Can also be used as a switch to activate toys
- Ideas for use can be found at <http://www.ablenetinc.com/portals/0/pdfs/101-ideas-bigsbs.pdf>

Big Points



- Record a single message on these.
- Press to hear the message played back.
- Come in 6 different colours.
- The clear plastic lid can be removed and a photo/picture/symbol inserted.


Doo-zy switch 000010



- Recordable switching communicator
- Manual, carrying case, adaptor and cable
- Doo-zy can be accessed directly through touch or by proximity
- Additional feedback can be given by the switch using vibration and/or lights
- Has multiple uses e.g. switch, recorder, infra red remote, toy controller, environmental control
- doo-zy.co.uk

Etran Frame

www.liberator.co.uk



- Serial no 0000119 (2 different layouts) and Stand
- Eye pointing is potentially a simple and useful way of communicating for people who cannot physically access other systems.
- The ETRAN frame can be used with photos, symbols or even objects for communication.
- Further advice on using an ETRAN frame can be found by going to <http://www.callscotland.org.uk/Common-Assets/quick-guides/QG-1408113630.pdf>

Switches

Switches can be used to help someone with physical difficulties access a variety of things within their environment such as toys, lights, music, food mixers.

Some items can be purchased already switch adapted but others would need battery adaptors/mains controller/click-on to facilitate this.

Assessment would need to be undertaken with your service user to decide which switch best meets their needs.

Switches vary greatly in size and access method. A range of these are available through the AAC kit in your area.





- Large Pal Pad
- Green RTS 270
- Activated using very light touch/movement

Pal Pad

- Pal Pad Small
- Red RTS 269
- Activated using very light touch/small movement



Specs switch

www.ablenetinc.com



RTS 265

RTS 266

screwdriver

strap

2 plastic discs

Buddy Button (Single switch)

www.ablenet.com


- Green RTS 272
- Yellow RTS 271
- 273
yellow/red/blue/green
- Sticky attachment rectangle



Microlight switch

www.ablenet.com

RTS 273
switch



Computer switch Interface

www.ablenet.com

1 usb cable and switch



This is a “plug-and-play” switch interface.

There is no software set up required

It can support up to 5 different functions using 5 switches.

Can be used with any switch accessible computer game/activity e.g. www.helpkidzlearn.com

Switch 2 Scan

www.pretorianuk.com



- Serial no. 000007, usb port connector and instruction leaflet
- Used to make an ipad switch accessible using scanning method.

Switch Latch and Timer



- This is great for a person using a single battery operated device
- There are 3 modes of control that determine how much time the device operates for once the switch is activated – minutes, seconds, on/off

Click-On



- This is a mains controller. It can operate two mains devices at the same time.
- It has different options available for operating the selected object/toy.
- More information from www.inclusive.co.uk

Makaton Resources



- Makaton is a language programme that uses signs and symbols to help people to communicate.
- These are some examples of resources available.
- More information and some free signing resources are available from www.makaton.org